



# «SMART» SYSTEMS:

a new generation of  
Thulium fiber lasers



*UROLAS+*

*UROLAS+  
PREMIUM*



# UROLAS<sup>+</sup>



## Thulium Fiber



## Lithotripsy

40 

## Special



«MRP» mode –



«Fine» dusting –



«Ultra» fragmentation –



«Dissect» enucleation mode –



«Bloodless» coagulation mode –



Tissue sensor mode –

## IRE SURGICAL FIBER

Design options:

• Single use • Multiple use

Available diameters

150  $\mu$ m



200  $\mu$ m




365  $\mu$ m



550  $\mu$ m



940  $\mu$ m



## One Push Connector

## laser «Smart»



**Lithotripsy**



**Soft tissues**

UROLAS<sup>+</sup>  
PREMIUM




## features

- pulse setting to minimize retropulsion
- ultra-fast fragmentation into micro-fragments
- breaking into large fragments for extraction
- thermo-mechanical dissection of tissues
- the most efficient coagulation mode
- tissue/stone detection




## Technical features

Standard  
network  
connection



Air-cooling



Regular  
maintenance  
is not required




4 times more compact  
and lighter than Ho:  
YAG high and medium  
power lasers



## Modulated pulses


Modulated pulse settings of **Urolase+** and **Urolase+ Premium** laser devices allow lithotripsy in different modes: from crushing «into dust» to breaking into large fragments for lithoextraction and lithoevacuation.

1




The new «**Fine**» **dusting** mode allows the surgeon to crush stones into fine dust at high speed.

2



The special «**Ultra**» **pulse** fragmentation mode instantly breaks down the densest stones into large fragments for subsequent lithoexcavation.

3




**MRP\*** **mode** - minimal stonedisplacementduring crushing, compared to holmium lasers and with standard pulses.



## Tissue Sensor – tissue/stone detection

**Tissue Sensor** is an innovative development of our company aimed at **absolute maximization of safety** during stone crushing.

This technology is designed to eliminate accidental exposure of soft tissues to laser radiation during lithotripsy.



The principle of the Tissue Sensor is that the laser detects which tissue (hard or soft) is in front of the surgical fiber tip.

Thus, during lithotripsy, the laser **automatically stops radiation** when it is pointed at the soft tissues, eliminating the risk of damage and perforation.



# Soft tissues




## Two types of enucleation in one device

The **Urolase+ Premium** laser device has two types of enucleation:

### 1 «Dissect» mode enucleation

- Adenomatous tissue dissection is the same as the HoLEP procedure
- Haemostatic properties are by far superior to those of HoLEP
- No carbonization




### 2 Classic thulium fiber enucleation – ThuFLEP

- Effective vaporization of soft tissues
- Precise work due to minimal depth of penetration
- No blood loss due to high level of hemostasis



## Coagulation Mode «Bloodless»

**Urolase+ Premium** has a unique pulse mode for coagulation. Due to its wide area of action, this mode allows effective coagulation of the postoperative area from a short distance.





# Accessories

## New wireless radiation activation pedal



2 in 1

In addition to the wireless connection, it is also possible to connect the pedal by wire, which is included in the kit.

## Urolase Cart\* Laser Trolley



## Technical characteristics

	UROLASE		UROLASE PREMIUM	
Wavelength, $\mu\text{m}$	1,94		1,94	
Laser type	Tm fiber		Tm fiber	
Operating mode	Pulsed	CW	Pulsed	CW
Maximum power, W	40		70	
Energy in pulse, J	0.02...6	-	0.02...6	-
Frequency, Hz	2000	-	3500	-
Cooling system	Air		Air	
Power supply voltage, V	220 $\pm$ 10 %		220 $\pm$ 10 %	
Network frequency, Hz	50...60		50...60	
Power consumption, V*A not more than	1600		1600	
Dimensions L*W*H, mm	606 x 526 x 314		606 x 526 x 314	
Weight, kg	45		45	

\*Urolase Cart is not included in the basic package of devices



# WORLD LEADER IN THE LASER INDUSTRY

---

VPG LaserOne LLC (formerly IRE-Polus LLC) is a vertically integrated company established by an outstanding Soviet scientist, Valentin Pavlovich Gapontsev, the founder of the international scientific and technical IPG Photonics Corporation.

VPG LaserOne is a globally recognized leader in the field of fiber lasers and amplifiers, as well as devices and systems based on them. Drawing on deep expertise and decades of experience in laser equipment production, VPG LaserOne LLC designs and supplies medical laser devices and surgical fiber for a wide range of applications.

VPG LaserOne develops advanced medical laser devices through a full-cycle process that includes device engineering, development of clinical application protocols, in-vitro research in its proprietary laboratory and clinical trials conducted in collaboration with leading clinical centers.



VPG LASERONE LLC  
[www.vpgmeds.com](http://www.vpgmeds.com)



+971 50 764 2603  
[sales@vpglaser.com](mailto:sales@vpglaser.com)



DATE OF  
ESTABLISHMENT  
**1992**

---



**15**  
CLINICAL CENTERS FOR  
IN-VITRO AND IN-VIVO  
STUDIES

---



**>1 million**  
PATIENTS TREATED WITH  
VPG LASERS IN 2024

---



**>2000**  
MEDICAL LASER SYSTEMS  
INSTALLED WORLDWIDE  
SINCE 2017